

**Product Name: TAOK1 Rabbit Monoclonal Antibody**  
**Catalog #: AMRe03222**



## Summary

<b>Production Name</b>	TAOK1 Rabbit Monoclonal Antibody
<b>Description</b>	Recombinant Rabbit Monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IP
<b>Reactivity</b>	Human,Mouse,Rat

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal Antibody
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
<b>Purification</b>	Affinity Purified

## Immunogen

<b>Gene Name</b>	TAOK1
<b>Alternative Names</b>	PSK2; TAO1; KFC-B; MARKK; PSK-2; hKFC-B; hTAOK1; MAP3K16
<b>Gene ID</b>	57551
<b>SwissProt ID</b>	Q7L7X3

## Application

<b>Dilution Ratio</b>	WB: 1/500-1/1000 IP: 1/20
<b>Molecular Weight</b>	Calculated MW: 116 kDa; Observed MW: 116 kDa

## Background

**Product Name: TAOK1 Rabbit Monoclonal Antibody**  
**Catalog #: AMRe03222**




Serine/threonine-protein kinase involved in various processes such as p38/MAPK14 stress-activated MAPK cascade, DNA damage response and regulation of cytoskeleton stability.


## Research Area

-

## Image Data



Western blot analysis of TAOK1 in 3T3, Raw264.7 lysates using TAOK1 antibody.



Western blot analysis of TAOK1 in C6, Ramos lysates using TAOK1 antibody.

## Note

For research use only.