

**Product Name: Ubiquitin Conjugation Factor E4 B Rabbit
Monoclonal Antibody
Catalog #: AMRe03250**

Summary

Production Name	Ubiquitin Conjugation Factor E4 B Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF
Reactivity	Human, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	UBE4B
Alternative Names	E4; UFD2; HDNB1; UBOX3; UFD2A
Gene ID	10277
SwissProt ID	O95155

Application

Dilution Ratio	WB: 1/500-1/1000 IF: 1/50-1/200
Molecular Weight	Calculated MW: 146 kDa; Observed MW: 146 kDa

**Product Name: Ubiquitin Conjugation Factor E4 B Rabbit
Monoclonal Antibody
Catalog #: AMRe03250**




Background

Ubiquitin-protein ligase that probably functions as an E3 ligase in conjunction with specific E1 and E2 ligases (By similarity).


Research Area

-


Image Data



Western blot analysis of UBE4B in PC-12, C6 lysates using Ubiquitin Conjugation Factor E4 B antibody.



Western blot analysis of UBE4B in Hela, A549, HL-60, U2OS, U897MG lysates using UBE4B antibody.



Immunocytochemistry analysis of UBE4B (green) in MCF-7 using UBE4B antibody, and DAPI(blue)



**Product Name: Ubiquitin Conjugation Factor E4 B Rabbit
Monoclonal Antibody
Catalog #: AMRe03250**

Note

For research use only.