

**Product Name: Phospho-PAK4/5/6  
(Ser474/Ser560/Ser602) Rabbit Monoclonal Antibody  
Catalog #: AMRe03279**

---

## Summary

<b>Production Name</b>	Phospho-PAK4/5/6 (Ser474/Ser560/Ser602) Rabbit Monoclonal Antibody
<b>Description</b>	Recombinant Rabbit Monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB
<b>Reactivity</b>	Human,Mouse,Rat

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Phosphorylated
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal Antibody
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
<b>Purification</b>	Affinity Purified

## Immunogen

<b>Gene Name</b>	PAK4/PAK5/PAK6
<b>Alternative Names</b>	PAK4; PAK5; PAK6
<b>Gene ID</b>	10298/57144/56924
<b>SwissProt ID</b>	O96013/Q9P286/Q9NQ5

## Application

<b>Dilution Ratio</b>	WB: 1/500-1/1000
<b>Molecular Weight</b>	Calculated MW: 64,80,75 kDa; Observed MW: 64,80,75 kDa

**Product Name: Phospho-PAK4/5/6  
(Ser474/Ser560/Ser602) Rabbit Monoclonal Antibody  
Catalog #: AMRe03279**

---


## Background

Serine/threonine protein kinase that plays a role in a variety of different signaling pathways including cytoskeleton regulation, cell migration, growth, proliferation or cell survival.

## Research Area

Signal Transduction

## Image Data



Western blot analysis of Phospho-PAK4 + PAK5 + PAK6 (Ser474 + Ser560 + Ser602) in HeLa, A549, HL-60, U2OS, U87-MG lysates using Phospho-PAK4/5/6 (Ser474/Ser560/Ser602) antibody.

## Note

For research use only.