

Product Name: Desmoglein 2 Rabbit Monoclonal Antibody
Catalog #: AMRe03983

Summary

Production Name	Desmoglein 2 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P
Reactivity	Human, Mouse, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
Purification	Affinity Purified

Immunogen

Gene Name	DSG2
Alternative Names	desmoglein-2; DSG2; HDGC; ARVD10; CDHF5; CMD1BB
Gene ID	1829
SwissProt ID	Q14126

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100
Molecular Weight	Calculated MW:122 kDa;Observed MW: 160 kDa

Product Name: Desmoglein 2 Rabbit Monoclonal Antibody
Catalog #: AMRe03983


Background

Component of intercellular desmosome junctions. Involved in the interaction of plaque proteins and intermediate filaments mediating cell-cell adhesion.

Research Area

Signal Transduction

Image Data



Western blot analysis of Desmoglein 2 in HepG2 lysates using Desmoglein 2 antibody.

Note

For research use only.