

Product Name: Plexin A1 Rabbit Monoclonal Antibody
Catalog #: AMRe04017



Summary

Production Name	Plexin A1 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,FC
Reactivity	Human, Mouse, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
Purification	Affinity Purified

Immunogen

Gene Name	PLXNA1
Alternative Names	NOV; NOV; Plexin-A1; Plexin1; PlexinA1; PLXN1; Plxna1
Gene ID	5361
SwissProt ID	Q9UIW2

Application

Dilution Ratio	WB: 1/500-1/1000 FC: 1/50-1/100
Molecular Weight	Calculated MW:211 kDa;Observed MW: 211 kDa

Background

Product Name: Plexin A1 Rabbit Monoclonal Antibody
Catalog #: AMRe04017




Plays a role in axon guidance, invasive growth and cell migration.

Research Area

Neuroscience

Image Data



Western blot analysis of Plexin A1 in K562 lysates using Plexin A1 antibody.

Note

For research use only.