

Product Name: Integrin alpha V Rabbit Monoclonal Antibody
Catalog #: AMRe04041



Summary

Production Name	Integrin alpha V Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF, IP
Reactivity	Human, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	ITGAV
Alternative Names	ITGAV; MSK8; VNRA; Integrin alpha-V; Vitronectin receptor subunit alpha; CD antigen CD51
Gene ID	3685
SwissProt ID	P06756

Application

Dilution Ratio	WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20
Molecular Weight	Calculated MW: 116 kDa; Observed MW: 130 kDa

Product Name: Integrin alpha V Rabbit Monoclonal Antibody
Catalog #: AMRe04041




Background

The alpha-V (ITGAV) integrins are receptors for vitronectin, cytotactin, fibronectin, fibrinogen, laminin, matrix metalloproteinase-2, osteopontin, osteomodulin, prothrombin, thrombospondin and vWF. They recognize the sequence R-G-D in a wide array of ligands. In case of HIV-1 infection, the interaction with extracellular viral Tat protein seems to enhance angiogenesis in Kaposi's sarcoma lesions.


Research Area

Cardiovascular

Image Data



Immunocytochemistry analysis of Integrin alpha V (green) in A549 using Integrin alpha V antibody, and DAPI(blue).



Western blot analysis of Integrin alphaV in C6, HeLa lysates using Integrin alphaV antibody.

Note

For research use only.