

Product Name: DAZAP1 Rabbit Monoclonal Antibody
Catalog #: AMRe04071



Summary

Production Name	DAZAP1 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	DAZAP1
Alternative Names	Deleted in azoospermia-associated protein 1
Gene ID	26528
SwissProt ID	Q96EP5

Application

Dilution Ratio	WB: 1/500-1/1000 IP: 1/20
Molecular Weight	Calculated MW: 43 kDa; Observed MW: 43 kDa

Background

Product Name: DAZAP1 Rabbit Monoclonal Antibody
Catalog #: AMRe04071




DAZ-associated protein 1

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of DAZAP1 in K562, C6, 3T3, HeLa lysates using DAZAP1 antibody.

Note

For research use only.