

Product Name: E2F4 Rabbit Monoclonal Antibody
Catalog #: AMRe04075



Summary

Production Name	E2F4 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	E2F4
Alternative Names	E2F4; Transcription factor E2F4; E2F-4
Gene ID	1874
SwissProt ID	Q16254

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 44 kDa; Observed MW: 65 kDa

Background

Product Name: E2F4 Rabbit Monoclonal Antibody
Catalog #: AMRe04075




Transcription activator that binds DNA cooperatively with DP proteins through the E2 recognition site, 5'-TTTC[CG]CGC-3' found in the promoter region of a number of genes whose products are involved in cell cycle regulation or in DNA replication.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of E2F4 in K562 lysates using E2F4 antibody.

Note

For research use only.