

**Product Name: Phospho-CREB (Ser133) (11Z4) Rabbit  
Monoclonal Antibody  
Catalog #: AMRe05881**

---



## Summary

|                        |                                                         |
|------------------------|---------------------------------------------------------|
| <b>Production Name</b> | Phospho-CREB (Ser133) (11Z4) Rabbit Monoclonal Antibody |
| <b>Description</b>     | Rabbit Monoclonal Antibody                              |
| <b>Host</b>            | Rabbit                                                  |
| <b>Application</b>     | WB                                                      |
| <b>Reactivity</b>      | Human,Mouse,Rat                                         |

## Performance

|                     |                                                                                                              |
|---------------------|--------------------------------------------------------------------------------------------------------------|
| <b>Conjugation</b>  | Unconjugated                                                                                                 |
| <b>Modification</b> | Phospho Antibody                                                                                             |
| <b>Isotype</b>      | IgG                                                                                                          |
| <b>Clonality</b>    | Monoclonal                                                                                                   |
| <b>Form</b>         | Liquid                                                                                                       |
| <b>Storage</b>      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                     |
| <b>Buffer</b>       | Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA. |
| <b>Purification</b> | Affinity purification                                                                                        |

## Immunogen

|                          |                                                                                              |
|--------------------------|----------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | CREB1                                                                                        |
| <b>Alternative Names</b> |                                                                                              |
| <b>Gene ID</b>           | 1385.0                                                                                       |
| <b>SwissProt ID</b>      | P16220.A synthetic phosphopeptide corresponding to residues surrounding Ser133 of human CREB |

## Application

|                         |                   |
|-------------------------|-------------------|
| <b>Dilution Ratio</b>   | WB: 1:1000-1:5000 |
| <b>Molecular Weight</b> | 37kDa             |

**Product Name: Phospho-CREB (Ser133) (11Z4) Rabbit  
Monoclonal Antibody  
Catalog #: AMRe05881**


---



## Background

## Research Area

## Image Data



Western blot detection of Phospho-CREB (Ser133) in CHO-K1 cell lysates using Phospho-CREB (Ser133) antibody(1:500 diluted).

## Note

For research use only.