

**Product Name: Phospho-Jun/JunD (Ser73/Ser100)
(12F13) Rabbit Monoclonal Antibody**
Catalog #: AMRe05932

Summary

Production Name	Phospho-Jun/JunD (Ser73/Ser100) (12F13) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Phospho Antibody
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	JUN/JUND
Alternative Names	
Gene ID	3725/3727
SwissProt ID	P05412/P17535.A synthetic phosphopeptide corresponding to residues surrounding Ser73 of human Jun

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	35,36kDa


**Product Name: Phospho-Jun/JunD (Ser73/Ser100)
(12F13) Rabbit Monoclonal Antibody**
Catalog #: AMRe05932



Background

Research Area

Image Data



Western blot detection of Jun (phospho-Ser73)/JunD (phospho-Ser100) in A549,HL-60,U251,U87-MG cell lysates using Jun (phospho-Ser73)/JunD (phospho-Ser100) antibody(1:1000 diluted).

Note

For research use only.