

Product Name: Phospho-Smad1 (Ser463/Ser465) (13F8)
Rabbit Monoclonal Antibody
Catalog #: AMRe06008

Summary

Production Name	Phospho-Smad1 (Ser463/Ser465) (13F8) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Phospho Antibody
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% New type preservative N and 0.05% BSA
Purification	Affinity purification

Immunogen

Gene Name	SMAD1
Alternative Names	
Gene ID	4086.0
SwissProt ID	Q15797.

Application


Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	52 kDa

Background

Product Name: Phospho-Smad1 (Ser463/Ser465) (13F8)
Rabbit Monoclonal Antibody
Catalog #: AMRe06008

Research Area

Image Data



Western blot analysis of Phospho-Smad1 (Ser463/Ser465) in A549, U2OS, C6 lysates using Phospho-Smad1 (Ser463/Ser465) antibody.

Note

For research use only.