

Product Name: Phospho-Smad3 (Ser423/425) (12O16)
Rabbit Monoclonal Antibody
Catalog #: AMRe06012



Summary

Production Name	Phospho-Smad3 (Ser423/425) (12O16) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Phospho Antibody
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	SMAD3
Alternative Names	
Gene ID	4088.0
SwissProt ID	P84022.A synthetic phosphopeptide corresponding to residues surrounding Ser423/Ser425 of human Smad3

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	48kDa


Product Name: Phospho-Smad3 (Ser423/425) (12016)
Rabbit Monoclonal Antibody
Catalog #: AMRe06012



Background

Research Area

Image Data



Western blot detection of Phospho-Smad3 (Ser423/425) in Jurkat,C6,CHO-K1,HeLa cell lysates using Phospho-Smad3 (Ser423/425) antibody(1:1000 diluted).

Note

For research use only.