

Product Name: Aconitase 2 (15R17) Rabbit Monoclonal Antibody
Catalog #: AMRe06512



Summary

Production Name	Aconitase 2 (15R17) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	ACO2
Alternative Names	ACO2; Aconitase2; Aconitate hydratase; ACONM; Citrate hydro lyase; ICRD;
Gene ID	50.0
SwissProt ID	Q99798.A synthetic peptide of human Aconitase 2

Application

Dilution Ratio	WB: 1:2000-1:20000
Molecular Weight	85kDa

Background


Product Name: Aconitase 2 (15R17) Rabbit Monoclonal Antibody
Catalog #: AMRe06512



Catalyzes the isomerization of citrate to isocitrate via cis-aconitate. Catalyzes the isomerization of citrate to isocitrate via cis- aconitate.

Research Area

Image Data



Western blot detection of Aconitase 2 in HL-60,C6 using Aconitase 2 antibody(1:1000 diluted)

Note

For research use only.