

Product Name: AMPK alpha (14Y14) Rabbit Monoclonal Antibody
Catalog #: AMRe06840

Summary

Production Name	AMPK alpha (14Y14) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	PRKAA2
Alternative Names	
Gene ID	5563.0
SwissProt ID	P54646.Recombinant protein of human AMPK alpha 2

Application


Dilution Ratio	WB: 1:1000
Molecular Weight	64,62kDa

Background

Product Name: AMPK alpha (14Y14) Rabbit Monoclonal Antibody
Catalog #: AMRe06840

Research Area

Image Data



Western blot detection of AMPK alpha 1/2 in Jurkat cell lysates using AMPK alpha 1/2 antibody(1:500 diluted).

Note

For research use only.