

Product Name: Amyloid Precursor Protein (5W2) Rabbit Monoclonal Antibody
Catalog #: AMRe06868

Summary

Production Name	Amyloid Precursor Protein (5W2) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	APP
Alternative Names	
Gene ID	351.0
SwissProt ID	P05067.A synthetic peptide of human Amyloid Precursor Protein

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	87kDa


Background

**Product Name: Amyloid Precursor Protein (5W2) Rabbit
Monoclonal Antibody
Catalog #: AMRe06868**



Research Area

Image Data



Western blot detection of Amyloid Precursor Protein in HeLa cell lysates using Amyloid Precursor Protein antibody(1:500 diluted).

Note

For research use only.