

Product Name: AP3M1 (16F4) Rabbit Monoclonal Antibody
Catalog #: AMRe06981

Summary

Production Name	AP3M1 (16F4) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	AP3M1
Alternative Names	
Gene ID	26985.0
SwissProt ID	Q9Y2T2.A synthetic peptide of human AP3M1

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	47kDa


Background

Product Name: AP3M1 (16F4) Rabbit Monoclonal Antibody
Catalog #: AMRe06981



Research Area

Image Data



Western blot detection of AP3M1 in K562,C6,3T3,HeLa cell lysates using AP3M1 antibody(1:1000 diluted).

Note

For research use only.