

Product Name: ATF4 (16F19) Rabbit Monoclonal Antibody
Catalog #: AMRe07275

Summary

Production Name	ATF4 (16F19) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol.
Purification	Affinity purification

Immunogen

Gene Name	ATF4
Alternative Names	
Gene ID	468.0
SwissProt ID	P18848.

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	39 kDa


Background

Product Name: ATF4 (16F19) Rabbit Monoclonal Antibody
Catalog #: AMRe07275



Research Area

Image Data



Western blot analysis of ATF4 in HeLa lysates using ATF4 antibody.

Note

For research use only.