

Product Name: Caspase 9 (15M16) Rabbit Monoclonal Antibody
Catalog #: AMRe07953



Summary

Production Name	Caspase 9 (15M16) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	Casp9
Alternative Names	
Gene ID	12371.0
SwissProt ID	Q8C3Q9.Recombinant protein of mouse Caspase-9

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	50kDa


Background

Product Name: Caspase 9 (15M16) Rabbit Monoclonal Antibody
Catalog #: AMRe07953



Research Area

Image Data



Western blot detection of Caspase-9 in 3T3 cell lysates using Caspase-9 antibody(1:1000 diluted).

Note

For research use only.