

Product Name: Caspase 9 (3P11) Rabbit Monoclonal Antibody
Catalog #: AMRe07954



Summary

Production Name	Caspase 9 (3P11) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	CASP9
Alternative Names	
Gene ID	842.0
SwissProt ID	P55211.Recombinant protein of human Caspase-9

Application


Dilution Ratio	WB: 1:1000
Molecular Weight	46kDa

Background

Product Name: Caspase 9 (3P11) Rabbit Monoclonal Antibody
Catalog #: AMRe07954

Research Area

Image Data



Western blot analysis of extracts from Mouse heart tissue using RM0070 at 1:1000.

Note

For research use only.