

Product Name: CBL (5W16) Rabbit Monoclonal Antibody
Catalog #: AMRe08037

Summary

Production Name	CBL (5W16) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	CBL
Alternative Names	
Gene ID	867.0
SwissProt ID	P22681.A synthetic peptide of human CBL

Application


Dilution Ratio	WB: 1:2000
Molecular Weight	100kDa

Background

Product Name: CBL (5W16) Rabbit Monoclonal Antibody
Catalog #: AMRe08037

Research Area

Image Data



Western blot detection of CBL in HeLa cell lysates using CBL antibody(1:1000 diluted).

Note

For research use only.