

Product Name: CDC123 (15X2) Rabbit Monoclonal Antibody
Catalog #: AMRe08496



Summary

Production Name	CDC123 (15X2) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	CDC123
Alternative Names	C10orf7; CDC123; Cell division cycle 123; D123; HT 1080; PZ32;
Gene ID	8872.0
SwissProt ID	O75794.A synthetic peptide of human CDC123

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	39kDa

Background


Product Name: CDC123 (15X2) Rabbit Monoclonal Antibody
Catalog #: AMRe08496



Required for S phase entry of the cell cycle. Required for S phase entry of the cell cycle.

Research Area

Image Data



Western blot detection of CDC123 in HeLa,A549,HL-60 cell lysates using CDC123 antibody(1:1000 diluted).

Note

For research use only.