

Product Name: CEP43 (7K4) Rabbit Monoclonal Antibody
Catalog #: AMRe08665

Summary

Production Name	CEP43 (7K4) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	CEP43
Alternative Names	FGFR1OP; FOP;
Gene ID	11116.0
SwissProt ID	O95684.A synthetic peptide of human FGFR1 Oncogene Partner

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	43kDa


Background

Product Name: CEP43 (7K4) Rabbit Monoclonal Antibody
Catalog #: AMRe08665

Required for anchoring microtubules to the centrosomes. Required for anchoring microtubules to the centrosomes (PubMed:16314388, PubMed:28659385). Required for ciliation (PubMed:28625565, PubMed:28659385).

Research Area

Image Data



Western blot detection of FGFR1 Oncogene Partner in A549,HL-60,U251,U87-MG cell lysates using FGFR1 Oncogene Partner antibody(1:1000 diluted).

Note

For research use only.