

Product Name: CHD4 (15O2) Rabbit Monoclonal Antibody
Catalog #: AMRe08746

Summary

Production Name	CHD4 (15O2) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	CHD4
Alternative Names	ATP dependent helicase CHD4; CHD4; Chromodomain helicase DNA binding protein 4; Mi 2 autoantigen 218 kDa protein; Mi 2b; Mi2 beta;
Gene ID	1108.0
SwissProt ID	Q14839.

Application

Dilution Ratio	WB 1:1000~1:5000
Molecular Weight	218kDa


Product Name: CHD4 (15O2) Rabbit Monoclonal Antibody
Catalog #: AMRe08746

Background

Component of the histone deacetylase NuRD complex which participates in the remodeling of chromatin by deacetylating histones. Component of the histone deacetylase NuRD complex which participates in the remodeling of chromatin by deacetylating histones.

Research Area

Image Data



Western blot analysis of CHD4 expression in 293T cell lysate.

Note

For research use only.