

Product Name: Collagen I alpha 1 (15S16) Rabbit Monoclonal Antibody
Catalog #: AMRe09208

Summary

Production Name	Collagen I alpha 1 (15S16) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	COL1A1
Alternative Names	
Gene ID	1277.0
SwissProt ID	P02452.Recombinant protein of human Collagen I

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	139kDa


Background

**Product Name: Collagen I alpha 1 (15S16) Rabbit
Monoclonal Antibody
Catalog #: AMRe09208**



Research Area

Image Data



Western blot detection of Collagen I in A549,HL-60,U2OS using Collagen I antibody(1:1000 diluted)

Note

For research use only.