

Product Name: CT A (6E19) Rabbit Monoclonal Antibody
Catalog #: AMRe09478

Summary

Production Name	CT A (6E19) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	PCYT1A
Alternative Names	
Gene ID	5130.0
SwissProt ID	P49585.A synthetic peptide of human PCYT1A

Application


Dilution Ratio	WB: 1:1000-1:5000
Molecular Weight	42kDa

Background

Product Name: CT A (6E19) Rabbit Monoclonal Antibody
Catalog #: AMRe09478

Research Area

Image Data



Western blot detection of CT A in MCF7,3T3,HeLa cell lysates using CT A antibody(1:1000 diluted).

Note

For research use only.