

Product Name: Cyclin T1 (2A7) Rabbit Monoclonal Antibody
Catalog #: AMRe09606

Summary

Production Name	Cyclin T1 (2A7) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	CCNT1
Alternative Names	
Gene ID	904.0
SwissProt ID	O60563.Recombinant protein of human Cyclin T1

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	81kDa


Background

Product Name: Cyclin T1 (2A7) Rabbit Monoclonal Antibody
Catalog #: AMRe09606



Research Area

Image Data



Western blot detection of Cyclin T1 in HeLa, Jurkat, K562 cell lysates using Cyclin T1 antibody (1:1000 diluted).

Note

For research use only.