

**Product Name: Cytokeratin 14 (12Q6) Rabbit
Monoclonal Antibody
Catalog #: AMRe09724**

Summary

Production Name	Cytokeratin 14 (12Q6) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol.
Purification	Affinity purification

Immunogen

Gene Name	KRT14
Alternative Names	
Gene ID	3861.0
SwissProt ID	P02533.

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	52 kDa


Background

**Product Name: Cytokeratin 14 (12Q6) Rabbit
Monoclonal Antibody
Catalog #: AMRe09724**



Research Area

Image Data



Western blot analysis of Cytokeratin 14 in Mouse skin lysates using Cytokeratin 14 antibody.

Note

For research use only.