

**Product Name: DAZAP1 (12Y15) Rabbit Monoclonal Antibody**  
**Catalog #: AMRe09798**

---

## Summary

|                        |                                           |
|------------------------|-------------------------------------------|
| <b>Production Name</b> | DAZAP1 (12Y15) Rabbit Monoclonal Antibody |
| <b>Description</b>     | Rabbit Monoclonal Antibody                |
| <b>Host</b>            | Rabbit                                    |
| <b>Application</b>     | WB                                        |
| <b>Reactivity</b>      | Human,Mouse,Rat                           |

## Performance

|                     |                                                                                                              |
|---------------------|--------------------------------------------------------------------------------------------------------------|
| <b>Conjugation</b>  | Unconjugated                                                                                                 |
| <b>Modification</b> | Unmodified                                                                                                   |
| <b>Isotype</b>      | IgG                                                                                                          |
| <b>Clonality</b>    | Monoclonal                                                                                                   |
| <b>Form</b>         | Liquid                                                                                                       |
| <b>Storage</b>      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                     |
| <b>Buffer</b>       | Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA. |
| <b>Purification</b> | Affinity purification                                                                                        |

## Immunogen

|                          |                                            |
|--------------------------|--------------------------------------------|
| <b>Gene Name</b>         | DAZAP1                                     |
| <b>Alternative Names</b> |                                            |
| <b>Gene ID</b>           | 26528.0                                    |
| <b>SwissProt ID</b>      | Q96EP5.Recombinant protein of human DAZAP1 |

## Application

|                         |                    |
|-------------------------|--------------------|
| <b>Dilution Ratio</b>   | WB: 1:2000-1:10000 |
| <b>Molecular Weight</b> | 43kDa              |

## Background

---


**Product Name: DAZAP1 (12Y15) Rabbit Monoclonal Antibody**  
**Catalog #: AMRe09798**

---



## Research Area

## Image Data



Western blot detection of DAZAP1 in K562,C6,3T3,HeLa cell lysates using DAZAP1 antibody(1:1000 diluted).

## Note

For research use only.