

Product Name: DAZAP1 (14U10) Rabbit Monoclonal Antibody
Catalog #: AMRe09799



Summary

Production Name	DAZAP1 (14U10) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	DAZAP1
Alternative Names	DAZ associated protein 1; Dazap1;
Gene ID	26528.0
SwissProt ID	Q96EP5.Recombinant protein of human DAZAP1

Application

Dilution Ratio	WB: 1:2000-1:10000
Molecular Weight	43kDa

Background


Product Name: DAZAP1 (14U10) Rabbit Monoclonal Antibody
Catalog #: AMRe09799



RNA-binding protein, which may be required during spermatogenesis. RNA-binding protein, which may be required during spermatogenesis.

Research Area

Image Data



Western blot detection of DAZAP1 in K562,C6,3T3,HeLa cell lysates using DAZAP1 antibody(1:1000 diluted).

Note

For research use only.