

Product Name: DNAJA4 (17J6) Rabbit Monoclonal Antibody
Catalog #: AMRe10068

Summary

Production Name	DNAJA4 (17J6) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	DNAJA4
Alternative Names	DNAJA4; MST104; MSTP104; PRO1472;
Gene ID	55466.0
SwissProt ID	Q8WW22.A synthetic peptide of human DNAJA4

Application

Dilution Ratio	WB: 1:1000-1:2000
Molecular Weight	45kDa

Background


Product Name: DNAJA4 (17J6) Rabbit Monoclonal Antibody
Catalog #: AMRe10068



Contains 1 CR-type zinc finger. Contains 1 J domain.

Research Area

Image Data



Western blot detection of DNAJA4 in C6,Rat Brain lysates using DNAJA4 antibody(1:500 diluted).

Note

For research use only.