

Product Name: FADD (18S3) Rabbit Monoclonal Antibody
Catalog #: AMRe10793

Summary

Production Name	FADD (18S3) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	
Alternative Names	FADD protein; GIG3; Growth inhibiting gene 3 protein; Mediator of receptor induced toxicity; MORT1;
Gene ID	
SwissProt ID	Q61160(mouse)(human:Q13158).A synthetic peptide of mouse FADD

Application

Dilution Ratio	WB: 1:1000-1:5000
Molecular Weight	23kDa


Product Name: FADD (18S3) Rabbit Monoclonal Antibody
Catalog #: AMRe10793

Background

Apoptotic adaptor molecule that recruits caspase-8 or caspase-10 to the activated Fas (CD95) or TNFR-1 receptors. The resulting aggregate called the death-inducing signaling complex (DISC) performs caspase-8 proteolytic activation. Active caspase-8 initiates the subsequent cascade of caspases mediating apoptosis.

Research Area

Image Data



Western blot detection of FADD in SHSY-5Y cell lysates using FADD antibody(1:1000 diluted).

Note

For research use only.