

Product Name: Fetuin A (6M11) Rabbit Monoclonal Antibody
Catalog #: AMRe10910



Summary

Production Name	Fetuin A (6M11) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	AHSG
Alternative Names	16G2; A2HS; AHS; Ahsg; FETUA; FETUB; FetuinA; FetuinB; Gugu; HSGA; IRL685; PRO2743;
Gene ID	197.0
SwissProt ID	P02765.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	39kDa

Product Name: Fetuin A (6M11) Rabbit Monoclonal Antibody
Catalog #: AMRe10910




Background

Promotes endocytosis, possesses opsonic properties and influences the mineral phase of bone. Shows affinity for calcium and barium ions. Promotes endocytosis, possesses opsonic properties and influences the mineral phase of bone. Shows affinity for calcium and barium ions.

Research Area

Image Data



Western blot analysis of Fetuin A expression in Human plasma lysate.

Note

For research use only.