

Product Name: Follistatin (5U17) Rabbit Monoclonal Antibody
Catalog #: AMRe11062



Summary

Production Name	Follistatin (5U17) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	FST
Alternative Names	Activin binding protein; Follistatin; FS; fst; Human follistatin gene;
Gene ID	10468.0
SwissProt ID	P19883.Recombinant protein of human Follistatin

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	38kDa

Background


Product Name: Follistatin (5U17) Rabbit Monoclonal Antibody
Catalog #: AMRe11062



Binds directly to activin and functions as an activin antagonist. Specific inhibitor of the biosynthesis and secretion of pituitary follicle stimulating hormone (FSH). Binds directly to activin and functions as an activin antagonist. Specific inhibitor of the biosynthesis and secretion of pituitary follicle stimulating hormone (FSH).

Research Area

Image Data



Western blot detection of Follistatin in Rat Brain lysates using Follistatin antibody(1:1000 diluted).

Note

For research use only.