

Product Name: FOXO3 (14Y16) Rabbit Monoclonal Antibody
Catalog #: AMRe11103

Summary

Production Name	FOXO3 (14Y16) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% New type preservative N and 0.05% BSA
Purification	Affinity purification

Immunogen

Gene Name	FOXO3
Alternative Names	
Gene ID	2309.0
SwissProt ID	O43524.

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	71 kDa


Background

Product Name: FOXO3 (14Y16) Rabbit Monoclonal Antibody
Catalog #: AMRe11103



Research Area

Image Data



Western blot analysis of FOXO3A in HeLa, A549 lysates using FOXO3 antibody.

Note

For research use only.