

Product Name: FUBP1 (9R5) Rabbit Monoclonal Antibody
Catalog #: AMRe11179

Summary

Production Name	FUBP1 (9R5) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	FUBP1
Alternative Names	
Gene ID	8880.0
SwissProt ID	Q96AE4.A synthetic peptide of human FUBP1

Application


Dilution Ratio	WB: 1:1000
Molecular Weight	68kDa

Background

Product Name: FUBP1 (9R5) Rabbit Monoclonal Antibody
Catalog #: AMRe11179

Research Area

Image Data



Western blot detection of FUBP1/FBP in K562,Rat Brain,HeLa cell lysates using FUBP1/FBP antibody(1:1000 diluted).

Note

For research use only.