

Product Name: GCLC (5X17) Rabbit Monoclonal Antibody
Catalog #: AMRe11352

Summary

Production Name	GCLC (5X17) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	GCLC
Alternative Names	GCS heavy chain; Gamma-ECS Gamma-glutamylcysteine synthetase; GLCL; GLCLC; GCLC
Gene ID	2729.0
SwissProt ID	P48506.A synthetic peptide of human GCLC

Application

Dilution Ratio	WB: 1:1000-1:5000
Molecular Weight	73kDa


Product Name: GCLC (5X17) Rabbit Monoclonal Antibody
Catalog #: AMRe11352

Background

Sulfur metabolism; glutathione biosynthesis; glutathione from L-cysteine and L-glutamate: step 1/2.

Research Area

Image Data



Western blot detection of GCLC in C6 cell lysates using GCLC antibody(1:1000 diluted).

Note

For research use only.