

Product Name: GCYA3 (12V4) Rabbit Monoclonal Antibody
Catalog #: AMRe11377

Summary

Production Name	GCYA3 (12V4) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	GUCY1A1
Alternative Names	GUCY1A3; GC-S-alpha-1; GUC1A3; GUCA3; GCS-alpha-3; GUCY1A1; GC-SA3; GCS-alpha-1; GUCSA3;
Gene ID	2982.0
SwissProt ID	Q02108.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	77kDa

Product Name: GCYA3 (12V4) Rabbit Monoclonal Antibody
Catalog #: AMRe11377




Background

Belongs to the adenylyl cyclase class-4/guanylyl cyclase family. Contains 1 guanylate cyclase domain.

Research Area

Image Data



Western blot analysis of GCYA3 expression in Human fetal kidney lysate.

Note

For research use only.