

**Product Name: Glucose 6 Phosphate Dehydrogenase
(5Y1) Rabbit Monoclonal Antibody**
Catalog #: AMRe11485



Summary

Production Name	Glucose 6 Phosphate Dehydrogenase (5Y1) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	G6PD
Alternative Names	
Gene ID	2539.0
SwissProt ID	P11413.A synthetic peptide of human Glucose 6 Phosphate Dehydrogenase

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	59kDa


Background

**Product Name: Glucose 6 Phosphate Dehydrogenase
(5Y1) Rabbit Monoclonal Antibody**
Catalog #: AMRe11485



Research Area

Image Data



Western blot detection of Glucose 6 Phosphate Dehydrogenase in HeLa,A549,U2OS using Glucose 6 Phosphate Dehydrogenase antibody(1:1000 diluted).

Note

For research use only.