

Product Name: GM130 (6E7) Rabbit Monoclonal Antibody
Catalog #: AMRe11527

Summary

Production Name	GM130 (6E7) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat,Monkey,Cow,Dog

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol.
Purification	Affinity purification

Immunogen

Gene Name	GOLGA2
Alternative Names	
Gene ID	2801.0
SwissProt ID	Q08379.

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	113 kDa


Background

Product Name: GM130 (6E7) Rabbit Monoclonal Antibody
Catalog #: AMRe11527



Research Area

Image Data



Western blot analysis of GM130 in HeLa, lysates using GM130 antibody.

Note

For research use only.