

Product Name: GRK3 (15U16) Rabbit Monoclonal Antibody
Catalog #: AMRe11781

Summary

Production Name	GRK3 (15U16) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	GRK3 {ECO:0000312 MIM:109636}
Alternative Names	ADRBK2; BARK2; GRK3;
Gene ID	157.0
SwissProt ID	P35626.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	80kDa


Product Name: GRK3 (15U16) Rabbit Monoclonal Antibody
Catalog #: AMRe11781

Background

Specifically phosphorylates the agonist-occupied form of the beta-adrenergic and closely related receptors. Specifically phosphorylates the agonist-occupied form of the beta-adrenergic and closely related receptors.

Research Area

Image Data



Western blot analysis of GRK3 expression in Ramos cell lysate.

Note

For research use only.