

Product Name: HNF 4 alpha (9Z3) Rabbit Monoclonal Antibody
Catalog #: AMRe12124



Summary

Production Name	HNF 4 alpha (9Z3) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% New type preservative N and 0.05% BSA
Purification	Affinity purification

Immunogen

Gene Name	Hnf4A
Alternative Names	
Gene ID	15378.0
SwissProt ID	P49698.

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	53 kDa


Background

Product Name: HNF 4 alpha (9Z3) Rabbit Monoclonal Antibody
Catalog #: AMRe12124



Research Area

Image Data



Western blot analysis of HNF4alpha in Raw264.7 lysates using HNF 4 alpha antibody.

Note

For research use only.