

Product Name: HOXA9 (4Y9) Rabbit Monoclonal Antibody
Catalog #: AMRe12170

Summary

Production Name	HOXA9 (4Y9) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	HOXA9
Alternative Names	HOX1; ABD-B; HOX1G; HOX1.7;
Gene ID	3205.0
SwissProt ID	P31269.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	30kDa


Product Name: HOXA9 (4Y9) Rabbit Monoclonal Antibody
Catalog #: AMRe12170

Background

Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis. Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis. Required for induction of E-selectin and VCAM-1, on the endothelial cells surface at sites of inflammation.

Research Area

Image Data



Western blot analysis of HOXA9 expression in 293T cell lysate.

Note

For research use only.