

Product Name: HPDL (3C9) Rabbit Monoclonal Antibody
Catalog #: AMRe12191



Summary

Production Name	HPDL (3C9) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	HPDL
Alternative Names	
Gene ID	84842.0
SwissProt ID	Q96IR7.A synthetic peptide of human HPDL

Application


Dilution Ratio	WB: 1:1000-1:5000
Molecular Weight	39kDa

Background

Product Name: HPDL (3C9) Rabbit Monoclonal Antibody
Catalog #: AMRe12191

Research Area

Image Data



Western blot detection of HPDL in HeLa,Caco2,A549 cell lysates using HPDL antibody(1:1000 diluted).

Note

For research use only.