

Product Name: Hsp70 1B (5Z4) Rabbit Monoclonal Antibody
Catalog #: AMRe12251

Summary

Production Name	Hsp70 1B (5Z4) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	Hspa1B
Alternative Names	
Gene ID	15511.0
SwissProt ID	P17879.Recombinant protein of mouse Hsp70

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	70kDa


Background

Product Name: Hsp70 1B (5Z4) Rabbit Monoclonal Antibody
Catalog #: AMRe12251



Research Area

Image Data



Western blot detection of Hsp70 in K562,Rat Brain,3T3,HeLa cell lysates using Hsp70 antibody(1:1000 diluted).

Note

For research use only.