

Product Name: Hsp90 (14P3) Rabbit Monoclonal Antibody
Catalog #: AMRe12259

Summary

Production Name	Hsp90 (14P3) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	
Alternative Names	HSP 90; HSP90A; HSP90B; LAP-2; HSPC1; HSPC2;
Gene ID	
SwissProt ID	P07900/P08238.A synthetic peptide of human Hsp90 alpha

Application

Dilution Ratio	WB: 1:1000-1:5000
Molecular Weight	85kDa


Background

Product Name: Hsp90 (14P3) Rabbit Monoclonal Antibody
Catalog #: AMRe12259

Molecular chaperone that promotes the maturation, structural maintenance and proper regulation of specific target proteins involved for instance in cell cycle control and signal transduction. Undergoes a functional cycle that is linked to its ATPase activity. This cycle probably induces conformational changes in the client proteins, thereby causing their activation.

Research Area

Image Data



Western blot detection of HSP90 in K562,C6,HeLa cell lysates using HSP90 antibody(1:1000 diluted).

Note

For research use only.