

Product Name: Iba1 (17Z18) Rabbit Monoclonal Antibody
Catalog #: AMRe12326

Summary

Production Name	Iba1 (17Z18) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	AIF1
Alternative Names	IBA1; IRT1; AIF-1; IRT-1;
Gene ID	199.0
SwissProt ID	P55008.

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	17kDa


Product Name: Iba1 (17Z18) Rabbit Monoclonal Antibody
Catalog #: AMRe12326

Background

Actin-binding protein that enhances membrane ruffling and RAC activation. Enhances the actin-bundling activity of LCP1. Binds calcium. Plays a role in RAC signaling and in phagocytosis. Actin-binding protein that enhances membrane ruffling and RAC activation. Enhances the actin-bundling activity of LCP1. Binds calcium. Plays a role in RAC signaling and in phagocytosis. May play a role in macrophage activation and function. Promotes the proliferation of vascular smooth muscle cells and of T-lymphocytes. Enhances lymphocyte migration. Plays a role in vascular inflammation.

Research Area

Image Data



Western blot analysis of Iba1 expression in U937 cell lysate.

Note

For research use only.