

Product Name: IL12 p40 (19Q5) Rabbit Monoclonal Antibody
Catalog #: AMRe12490



Summary

Production Name	IL12 p40 (19Q5) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	IL12B
Alternative Names	CLMF; CLMF p40; CLMF2; IL12B; NKSF; NKSF2;
Gene ID	3593.0
SwissProt ID	P29460.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	37kDa

Product Name: IL12 p40 (19Q5) Rabbit Monoclonal Antibody
Catalog #: AMRe12490




Background

Cytokine that can act as a growth factor for activated T and NK cells, enhance the lytic activity of NK/lymphokine-activated killer cells, and stimulate the production of IFN-gamma by resting PBMC. Cytokine that can act as a growth factor for activated T and NK cells, enhance the lytic activity of NK/lymphokine-activated killer cells, and stimulate the production of IFN-gamma by resting PBMC.

Research Area

Image Data



Western blot analysis of IL12 p40 expression in Human IL12 p40 recombinant protein.

Note

For research use only.