

Product Name: JNK (19C18) Rabbit Monoclonal Antibody
Catalog #: AMRe12841

Summary

Production Name	JNK (19C18) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	MAPK8
Alternative Names	
Gene ID	5599.0
SwissProt ID	P45983.Recombinant protein of human JNK1

Application


Dilution Ratio	WB: 1:1000
Molecular Weight	48,53kDa

Background

Product Name: JNK (19C18) Rabbit Monoclonal Antibody
Catalog #: AMRe12841

Research Area

Image Data



Western blot detection of JNK1/2/3 in Jurkat,Rat Brain,C6,CHO-K1,HeLa cell lysates using JNK1/2/3 antibody(1:1000 diluted).

Note

For research use only.