

Product Name: LDHB (13O6) Rabbit Monoclonal Antibody
Catalog #: AMRe13264

Summary

Production Name	LDHB (13O6) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	LDHB
Alternative Names	HEL S 281; LDH B; LDH H; LDH-B; LDH-H; LDHB; LDHBD; LDHH; TRG5;
Gene ID	3945.0
SwissProt ID	P07195.

Application

Dilution Ratio	WB 1:1000-1:5000
Molecular Weight	37kDa


Product Name: LDHB (13O6) Rabbit Monoclonal Antibody
Catalog #: AMRe13264

Background

Fermentation; pyruvate fermentation to lactate; (S)-lactate from pyruvate: step 1/1.

Research Area

Image Data



Western blot analysis of LDHB expression in HeLa cell lysate.

Note

For research use only.