

Product Name: LGALS3BP (7X1) Rabbit Monoclonal Antibody
Catalog #: AMRe13287



Summary

Production Name	LGALS3BP (7X1) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	LGALS3BP
Alternative Names	
Gene ID	3959.0
SwissProt ID	Q08380.Recombinant protein of human LGALS3BP

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	65kDa


Background

Product Name: LGALS3BP (7X1) Rabbit Monoclonal Antibody
Catalog #: AMRe13287



Research Area

Image Data



Western blot detection of LGALS3BP in HT1080, HepG2, A549 cell lysates using LGALS3BP antibody(1:1000 diluted).

Note

For research use only.